



# Prolonged Culturing of Human Hepatocytes in Human Plasma for P450 Induction and *In Vitro* Hepatotoxicity Studies

Qian Yang, Carol J. Loretz, Kirsten Amaral, Doshi Utkarsh, and Albert P. Li, *In Vitro* ADMET Laboratories Inc., Columbia, MD and Malden, MA

## Introduction

- *In vitro* toxicology studies are routinely performed with the target cells cultured in protein free medium
- *In vivo* target cells are exposed to toxicants in the presence whole plasma
- An ideal *in vitro* model therefore should consist of the target cells cultured in whole plasma
- We hypothesize that evaluation of drug toxicity with a physiologically relevant target cell cultured in human plasma may allow direct extrapolation of *in vitro* findings to human *in vivo*
- Based on this hypothesis, we initiated the evaluation of the ability of human plasma to support the culturing of human hepatocytes for *in vitro* hepatotoxicity studies

## Materials & Methods


**Cell Culture:** Hepatocytes from donor HH1062 (IVAL, Columbia, MD) were plated in both human plasma (minimally-modified human plasma (HPZ-A), IVAL) and protein free hepatocyte media (HIM, IVAL) for a total duration of 29 days. The hepatocytes were plated in 96-well collagen-coated tissue culture plates using Universal Cryopreservation Plating Medium (UCPM, IVAL). After 4 h of attachment, the medium was changed to HPZ-A and HIM containing 0.25 mg/mL Matrigel. Medium was changed to HPZ-A without Matrigel on the next day, and every 2 – 3 days afterwards for a culture duration of 29 days.

For CYP induction comparisons, 4 lots of human hepatocytes (HH1007, HH1026, HH1051, HH1053) were used. The hepatocytes were plated in 96-well plates as described above in both HIM and HPZ-A, with treatment of prototypical inducers beginning on day 5 after plating. The treatment duration was 3 days. Gene expression by RT-PCR was used to measure induction of CYP1A2 (by omeprazole), CYP2B6 (by phenobarbital), and CYP3A4 (by rifampin).

For long-term cytotoxicity, lot HH1062 was cultured in HPZ-A as described above, and treated with seven concentrations of aflatoxin B1 (concentrations ranging from 1 uM to 90 uM) on day 2 and day 10. Viability was quantified based on cellular ATP contents (ATPLite, Perkin-Elmer). Results are expressed as Relative Viability using the following equation:

$$\text{Relative Viability (\%)} = \frac{\text{ATP (Treatment)}}{\text{ATP (Solvent Control)}} \times 100$$

## Morphology of Human Hepatocytes Cultured in Human Plasma




## P450 Induction in Human Plasma

Fold induction in four lots of human hepatocytes cultured in protein free hepatocyte medium (HIM) and modified human plasma (HPZ-A) for CYP1A2, CYP2B6, and CYP3A4

Lots	Fold Induction (mRNA)					
	CYP1A2		CYP2B6		CYP3A4	
	HIM	HPZ-A	HIM	HPZ-A	HIM	HPZ-A
HH1007	3	6	3	2	9	20
HH1026	4	7	4	2	13	25
HH1051	4	5	7	4	10	18
HH1053	4	8	3	2	9	16

## In Vitro Hepatotoxicity Evaluation in Human Plasma

Application of Human Hepatocyte/Human Plasma Culture for the Evaluation of Acute and Subchronic Hepatotoxicity: Proof-of-Principle Study with Aflatoxin B1



## Summary and Conclusions

- We have successfully developed a 100% human plasma medium, HPZ-A with which human hepatocytes could be cultured for a prolonged period of >29 days.
- Human hepatocytes cultured in HPZ-A are responsive to CYPs 1A2, 2B6, and 3A4 induction, demonstrating the retention of intact P450 induction pathways
- Aflatoxin B1, which requires metabolic activation for hepatotoxicity, was found to be cytotoxic to human hepatocytes cultured in human plasma. Aflatoxin IC<sub>50</sub> value after a 10-day treatment duration was 19X lower than that observed after a 2-day treatment duration. **The results suggest that human hepatocytes cultured in human plasma could activate a protoxicant and can be used for the evaluation of chronic hepatotoxicity.**

Contact: Albert P. Li, Ph. D.  
lialbert@invitroadmet.com